

## Mouse Asgr1/Asialoglycoprotein receptor 1 ELISA Kit

Catalog No: #EK19519

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Asgr1/Asialoglycoprotein receptor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Asialoglycoprotein receptor 1
Target Name	Asialoglycoprotein receptor 1
Other Names	ASGP-R 1,Hepatic lectin 1,HL-1,Asgr-1
Accession No.	P34927
Uniprot	P34927
GeneID	11889;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.099ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Asgr1

Note: This product is for in vitro research use only