

Mouse Comp/Cartilage oligomeric matrix protein ELISA Kit



Catalog No: #EK19521

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Comp/Cartilage oligomeric matrix protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cartilage oligomeric matrix protein
Target Name	Cartilage oligomeric matrix protein
Other Names	COMP
Accession No.	Q9R0G6
Uniprot	Q9R0G6
GeneID	12845;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.089ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Comp

Note: This product is for in vitro research use only