

## Mouse Csk/Tyrosine-protein kinase CSK ELISA Kit

Catalog No: #EK19522

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Csk/Tyrosine-protein kinase CSK ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Tyrosine-protein kinase CSK
Target Name	Tyrosine-protein kinase CSK
Other Names	C-Src kinase,Protein-tyrosine kinase MPK-2,p50CSK
Accession No.	P41241
Uniprot	P41241
GeneID	12988;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Csk

Note: This product is for in vitro research use only