

Rat Csk/Tyrosine-protein kinase CSK ELISA Kit

Catalog No: #EK19523

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Csk/Tyrosine-protein kinase CSK ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Tyrosine-protein kinase CSK
Target Name	Tyrosine-protein kinase CSK
Other Names	C-Src kinase
Accession No.	P32577
Uniprot	P32577
GeneID	315707;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.083ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Csk

Note: This product is for in vitro research use only