

Human CAST/Calpastatin ELISA Kit

Catalog No: #EK19524

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CAST/Calpastatin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Calpastatin
Target Name	Calpastatin
Other Names	Calpain inhibitor, Sperm BS-17 component
Accession No.	P20810
Uniprot	P20810
GeneID	831;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 30pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CAST

Note: This product is for in vitro research use only