

Rat Cast/Calpastatin ELISA Kit

Catalog No: #EK19526

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Rat Cast/Calpastatin ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Rt |
| Specificity | Natural and recombinant Rat Calpastatin |
| Target Name | Calpastatin |
| Other Names | Calpain inhibitor |
| Accession No. | P27321 |
| Uniprot | P27321 |
| GeneID | 25403; |

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:35pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cast

Note: This product is for in vitro research use only