

## Bovine CHGA/Chromogranin-A ELISA Kit

Catalog No: #EK19527

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CHGA/Chromogranin-A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Chromogranin-A
Target Name	Chromogranin-A
Other Names	CgA, Pituitary secretory protein I, SP-I
Accession No.	P05059
Uniprot	P05059
GeneID	281070;

## Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.78pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CHGA

Note: This product is for in vitro research use only