

Pig CHGA/Chromogranin-A ELISA Kit

Catalog No: #EK19528

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig CHGA/Chromogranin-A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Chromogranin-A
Target Name	Chromogranin-A
Other Names	CgA
Accession No.	P04404
Uniprot	P04404
GeneID	397540;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CHGA

Note: This product is for in vitro research use only