

# Chicken CALCA/Calcitonin gene-related peptide ELISA Kit

Catalog No: #EK19529

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken CALCA/Calcitonin gene-related peptide ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Calcitonin gene-related peptide
Target Name	Calcitonin gene-related peptide
Other Names	CGRP,CALC
Accession No.	P10286
Uniprot	P10286

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.063ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CALCA

Note: This product is for in vitro research use only