

Mouse Calca/Calcitonin gene-related peptide 1 ELISA Kit



Catalog No: #EK19530

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Calca/Calcitonin gene-related peptide 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Calcitonin gene-related peptide 1
Target Name	Calcitonin gene-related peptide 1
Other Names	Alpha-type CGRP, Calcitonin gene-related peptide I, CGRP-I, Calc
Accession No.	Q99JA0
Uniprot	Q99JA0
GeneID	12310;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.082ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Calca

Note: This product is for in vitro research use only