

## Mouse Adora1/Adenosine receptor A1 ELISA Kit

Catalog No: #EK19531

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Adora1/Adenosine receptor A1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Adenosine receptor A1
Target Name	Adenosine receptor A1
Accession No.	Q60612
Uniprot	Q60612
GeneID	11539;

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.059ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Adora1

Note: This product is for in vitro research use only