

Mouse Ctla4/Cytotoxic T-lymphocyte protein 4 ELISA Kit

Catalog No: #EK19533

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ctla4/Cytotoxic T-lymphocyte protein 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cytotoxic T-lymphocyte protein 4
Target Name	Cytotoxic T-lymphocyte protein 4
Other Names	Cytotoxic T-lymphocyte-associated antigen 4,CTLA-4,CD152,Cd152
Accession No.	P09793
Uniprot	P09793
GeneID	12477;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.084ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ctla4

Note: This product is for in vitro research use only