

## Bovine CDK4/Cyclin-dependent kinase 4 ELISA Kit

Catalog No: #EK19534

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CDK4/Cyclin-dependent kinase 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cyclin-dependent kinase 4
Target Name	Cyclin-dependent kinase 4
Other Names	Cell division protein kinase 4
Accession No.	Q32KY4
Uniprot	Q32KY4
GeneID	510618;

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.058ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CDK4

Note: This product is for in vitro research use only