

## Mouse Cdk4/Cyclin-dependent kinase 4 ELISA Kit

Catalog No: #EK19535

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cdk4/Cyclin-dependent kinase 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cyclin-dependent kinase 4
Target Name	Cyclin-dependent kinase 4
Other Names	CRK3, Cell division protein kinase 4, PSK-J3, Crk3
Accession No.	P30285
Uniprot	P30285
GeneID	12567;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.3pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cdk4

Note: This product is for in vitro research use only