

Rat Cdk4/Cyclin-dependent kinase 4 ELISA Kit

Catalog No: #EK19536

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cdk4/Cyclin-dependent kinase 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cyclin-dependent kinase 4
Target Name	Cyclin-dependent kinase 4
Other Names	Cell division protein kinase 4,PSK-J3
Accession No.	P35426
Uniprot	P35426
GeneID	94201;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.4pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cdk4

Note: This product is for in vitro research use only