

Bovine BPI/Bactericidal permeability-increasing protein ELISA Kit

Catalog No: #EK19537

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine BPI/Bactericidal permeability-increasing protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Bactericidal permeability-increasing protein
Target Name	Bactericidal permeability-increasing protein
Other Names	BPI
Accession No.	P17453
Uniprot	P17453
GeneID	280734;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.432ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: BPI

Note: This product is for in vitro research use only