

Mouse Bpi/Bactericidal permeability-increasing protein ELISA Kit



Catalog No: #EK19538

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Bpi/Bactericidal permeability-increasing protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bactericidal permeability-increasing protein
Target Name	Bactericidal permeability-increasing protein
Other Names	BPI
Accession No.	Q67E05
Uniprot	Q67E05
GeneID	329547;

Application Details

Assay Type: Sandwich
Detect Range: 31.2-2000pg/mL
Sensitivity: 15.82pg/mL
Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul
Assay Time: 3-5h
Detection wavelength: 450 nm
gene_full_name: Bpi

Note: This product is for in vitro research use only