

# Rat Bpi/Bactericidal permeability-increasing protein ELISA Kit



Catalog No: #EK19539

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Bpi/Bactericidal permeability-increasing protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Bactericidal permeability-increasing protein
Target Name	Bactericidal permeability-increasing protein
Other Names	BPI
Accession No.	Q6AXU0
Uniprot	Q6AXU0
GeneID	296321;

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.087ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Bpi

Note: This product is for in vitro research use only