

Mouse Ap1m1/AP-1 complex subunit mu-1 ELISA Kit

Catalog No: #EK19540

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ap1m1/AP-1 complex subunit mu-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse AP-1 complex subunit mu-1
Target Name	AP-1 complex subunit mu-1
Other Names	AP-mu chain family member mu1A,Adaptor protein complex AP-1 subunit mu-1,Adaptor-related protein complex 1 subunit mu-1,Clathrin assembly protein complex 1 mu-1 medium chain 1,Clathrin coat assembly protein AP47,Clathrin coat-associated protein AP47,Golgi adaptor HA1/AP1 adaptin mu-1 subunit,Mu-adaptin 1,Mu1A-adaptin,Cltnm
Accession No.	P35585
Uniprot	P35585
GeneID	11767;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.053ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ap1m1

Note: This product is for in vitro research use only