

Mouse Cotl1/Coactosin-like protein ELISA Kit

Catalog No: #EK19541

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cotl1/Coactosin-like protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Coactosin-like protein
Target Name	Coactosin-like protein
Other Names	Clp
Accession No.	Q9CQI6
Uniprot	Q9CQI6
GeneID	72042;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cotl1

Note: This product is for in vitro research use only