

Mouse Asph/Aspartyl/asparaginyl beta-hydroxylase ELISA Kit



Catalog No: #EK19544

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Asph/Aspartyl/asparaginyl beta-hydroxylase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Aspartyl/asparaginyl beta-hydroxylase
Target Name	Aspartyl/asparaginyl beta-hydroxylase
Other Names	Aspartate beta-hydroxylase,ASP beta-hydroxylase,Peptide-aspartate beta-dioxygenase,Bah
Accession No.	Q8BSY0
Uniprot	Q8BSY0
GeneID	65973;

Application Details

Assay Type: Sandwich Detect Range:31.2-2000pg/mL Sensitivity:15.6pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Asph

Note: This product is for in vitro research use only