

Mouse Abca1/ATP-binding cassette sub-family A member 1 ELISA Kit

Catalog No: #EK19545

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Abca1/ATP-binding cassette sub-family A member 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse ATP-binding cassette sub-family A member 1
Target Name	ATP-binding cassette sub-family A member 1
Other Names	ATP-binding cassette transporter 1,ABC-1,Abc1
Accession No.	P41233
Uniprot	P41233
GeneID	11303;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.069ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Abca1

Note: This product is for in vitro research use only