

Mouse Aadat/Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial ELISA Kit

Catalog No: #EK19553

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Aadat/Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial
Target Name	Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial
Other Names	KAT/AadAT,2-aminoadipate aminotransferase,2-aminoadipate transaminase,Alpha-aminoadipate aminotransferase,AadAT,Kynurenine aminotransferase II,Kynurenine--oxoglutarate aminotransferase II,Kynurenine--oxoglutarate transaminase 2,Kynurenine--oxoglutarate transaminase II,Kat2
Accession No.	Q9WVM8
Uniprot	Q9WVM8
GeneID	23923;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
 Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aadat

Note: This product is for in vitro research use only