

# Rat Atp1a1/Sodium/potassium-transporting ATPase subunit alpha-1 ELISA Kit



Catalog No: #EK19554

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Atp1a1/Sodium/potassium-transporting ATPase subunit alpha-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Sodium/potassium-transporting ATPase subunit alpha-1
Target Name	Sodium/potassium-transporting ATPase subunit alpha-1
Other Names	Na(+)/K(+) ATPase alpha-1 subunit,Sodium pump subunit alpha-1
Accession No.	P06685
Uniprot	P06685
GeneID	24211;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.083ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Atp1a1

Note: This product is for in vitro research use only