

Rat Atp2b1/Plasma membrane calcium-transporting ATPase 1 ELISA Kit

Catalog No: #EK19555

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Atp2b1/Plasma membrane calcium-transporting ATPase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Plasma membrane calcium-transporting ATPase 1
Target Name	Plasma membrane calcium-transporting ATPase 1
Other Names	PMCA1,Plasma membrane calcium ATPase isoform 1,Plasma membrane calcium pump isoform 1
Accession No.	P11505
Uniprot	P11505
GeneID	29598;

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:28pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Atp2b1

Note: This product is for in vitro research use only