

Human ATP2B2/Plasma membrane calcium-transporting ATPase 2 ELISA Kit

Catalog No: #EK19556

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ATP2B2/Plasma membrane calcium-transporting ATPase 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Plasma membrane calcium-transporting ATPase 2
Target Name	Plasma membrane calcium-transporting ATPase 2
Other Names	PMCA2, Plasma membrane calcium ATPase isoform 2, Plasma membrane calcium pump isoform 2, PMCA2
Accession No.	Q01814
Uniprot	Q01814
GeneID	491;

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.6pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ATP2B2

Note: This product is for in vitro research use only