

Mouse Atp2b2/Plasma membrane calcium-transporting ATPase 2 ELISA Kit

Catalog No: #EK19557

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse Atp2b2/Plasma membrane calcium-transporting ATPase 2 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Plasma membrane calcium-transporting ATPase 2 |
| Target Name | Plasma membrane calcium-transporting ATPase 2 |
| Other Names | PMCA2, Plasma membrane calcium ATPase isoform 2, Plasma membrane calcium pump isoform 2, Pmca2 |
| Accession No. | Q9R0K7 |
| Uniprot | Q9R0K7 |
| GeneID | 11941; |

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.085ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Atp2b2

Note: This product is for in vitro research use only