

Rat Atp2b2/Plasma membrane calcium-transporting ATPase 2 ELISA Kit



Catalog No: #EK19558

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Atp2b2/Plasma membrane calcium-transporting ATPase 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Plasma membrane calcium-transporting ATPase 2
Target Name	Plasma membrane calcium-transporting ATPase 2
Other Names	PMCA2, Plasma membrane calcium ATPase isoform 2, Plasma membrane calcium pump isoform 2
Accession No.	P11506
Uniprot	P11506
GeneID	24215;

Application Details

Assay Type: Sandwich Detect Range: 0.132-20ng/mL Sensitivity: 0.061ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Atp2b2

Note: This product is for in vitro research use only