

Mouse Atp5b/ATP synthase subunit beta, mitochondrial ELISA Kit

Catalog No: #EK19559

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Atp5b/ATP synthase subunit beta, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse ATP synthase subunit beta, mitochondrial
Target Name	ATP synthase subunit beta, mitochondrial
Accession No.	P56480
Uniprot	P56480
GeneID	11947;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.161ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Atp5b

Note: This product is for in vitro research use only