

Mouse Atp5j/ATP synthase-coupling factor 6, mitochondrial ELISA Kit

Catalog No: #EK19560

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Atp5j/ATP synthase-coupling factor 6, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse ATP synthase-coupling factor 6, mitochondrial
Target Name	ATP synthase-coupling factor 6, mitochondrial
Other Names	ATPase subunit F6
Accession No.	P97450
Uniprot	P97450
GeneID	11957;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.43ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Atp5j

Note: This product is for in vitro research use only