

Mouse Cdx2/Homeobox protein CDX-2 ELISA Kit

Catalog No: #EK19562

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cdx2/Homeobox protein CDX-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Homeobox protein CDX-2
Target Name	Homeobox protein CDX-2
Other Names	Caudal-type homeobox protein 2, Cdx-2
Accession No.	P43241
Uniprot	P43241
GeneID	12591;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 40.5pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdx2

Note: This product is for in vitro research use only