

Mouse Cdh7/Cadherin-7 ELISA Kit

Catalog No: #EK19563

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cdh7/Cadherin-7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cadherin-7
Target Name	Cadherin-7
Accession No.	Q8BM92
Uniprot	Q8BM92
GeneID	241201;

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 8.08pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdh7

Note: This product is for in vitro research use only