

## Mouse Btnl2/Butyrophilin-like protein 2 ELISA Kit

Catalog No: #EK19565

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Btnl2/Butyrophilin-like protein 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Butyrophilin-like protein 2
Target Name	Butyrophilin-like protein 2
Other Names	Gm315,Ng9
Accession No.	O70355
Uniprot	O70355
GeneID	547431;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.08ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Btnl2

Note: This product is for in vitro research use only