

Rat Cd46/Membrane cofactor protein ELISA Kit

Catalog No: #EK19567

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cd46/Membrane cofactor protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Membrane cofactor protein
Target Name	Membrane cofactor protein
Other Names	CD46,Mcp
Accession No.	Q9Z0M4
Uniprot	Q9Z0M4
GeneID	29333;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Cd46

Note: This product is for in vitro research use only