

# Mouse Cirbp/Cold-inducible RNA-binding protein ELISA Kit



Catalog No: #EK19568

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Mouse Cirbp/Cold-inducible RNA-binding protein ELISA Kit         |
| Brief Description  | ELISA Kits   |
| Applications       | ELISA  |
| Species Reactivity | Ms   |
| Specificity        | Natural and recombinant Mouse Cold-inducible RNA-binding protein |
| Target Name        | Cold-inducible RNA-binding protein                               |
| Other Names        | A18 hnRNP, Glycine-rich RNA-binding protein CIRP, Cirp           |
| Accession No.      | P60824   |
| Uniprot            | P60824   |
| GeneID             | 12696;   |

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.086ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cirbp

Note: This product is for in vitro research use only