

Bovine ADRA1A/Alpha-1A adrenergic receptor ELISA Kit



Catalog No: #EK19569

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine ADRA1A/Alpha-1A adrenergic receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Alpha-1A adrenergic receptor
Target Name	Alpha-1A adrenergic receptor
Other Names	Alpha-1A adrenoreceptor,Alpha-1A adrenoceptor,Alpha-1C adrenergic receptor,ADRA1C
Accession No.	P18130
Uniprot	P18130
GeneID	282134;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.087ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ADRA1A

Note: This product is for in vitro research use only