

Human ADRA1A/Alpha-1A adrenergic receptor ELISA Kit



Catalog No: #EK19570

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ADRA1A/Alpha-1A adrenergic receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Alpha-1A adrenergic receptor
Target Name	Alpha-1A adrenergic receptor
Other Names	Alpha-1A adrenoreceptor,Alpha-1A adrenoceptor,Alpha-1C adrenergic receptor,Alpha-adrenergic receptor 1c,ADRA1C
Accession No.	P35348
Uniprot	P35348
GeneID	148;

Application Details

Assay Type: SandwichDetect Range:0.031-2ng/ml Sensitivity:0.0155ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:ADRA1A

Note: This product is for in vitro research use only