

Mouse Adra1a/Alpha-1A adrenergic receptor ELISA Kit

Catalog No: #EK19571

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Adra1a/Alpha-1A adrenergic receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Alpha-1A adrenergic receptor
Target Name	Alpha-1A adrenergic receptor
Other Names	Alpha-1A adrenoreceptor,Alpha-1A adrenoceptor,Alpha-1C adrenergic receptor,Adra1c
Accession No.	P97718
Uniprot	P97718
GeneID	11549;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.091ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Adra1a

Note: This product is for in vitro research use only