

Rat Adra1a/Alpha-1A adrenergic receptor ELISA Kit

Catalog No: #EK19572

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Adra1a/Alpha-1A adrenergic receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Alpha-1A adrenergic receptor
Target Name	Alpha-1A adrenergic receptor
Other Names	Alpha-1A adrenoreceptor, Alpha-1A adrenoceptor, Alpha-1C adrenergic receptor, Adra1c
Accession No.	P43140
Uniprot	P43140
GeneID	29412;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.184ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adra1a

Note: This product is for in vitro research use only