

GNA12 Antibody

Catalog No: #44665

Package Size: #44665-1 50ul #44665-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | GNA12 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Immunogen affinity purified |
| Applications | WB |
| Species Reactivity | Hu Ms Rt |
| Specificity | GNA12antibody detects endogenous levels of total?GNA12 |
| Immunogen Type | Peptide |
| Immunogen Description | A synthesized peptide derived from human GNA12 |
| Target Name | GNA12 |
| Other Names | G alpha 12 antibody G alpha-12 antibody G-protein subunit alpha-12 antibody gep antibody Gna12 antibody GNA12_HUMAN antibody guanine nucleotide binding protein (G protein) alpha 12 antibody Guanine nucleotide binding protein alpha 12 subunit antibody Guanine nucleotide-binding protein subunit alpha- |
| Accession No. | Swiss-Prot#: Q03113NCBI Gene ID:2768 |
| Uniprot | Q03113 |
| GeneID | 2768; |
| Calculated MW | 44 kDa |
| Concentration | 1.0mg/mL |
| Formulation | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Storage | Store at -20°C |

Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only