

Mouse Agrn/Agrin ELISA Kit

Catalog No: #EK19574

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Agrn/Agrin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Agrin
Target Name	Agrin
Other Names	Agrin
Accession No.	A2ASQ1
Uniprot	A2ASQ1
GeneID	11603;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Agrn

Note: This product is for in vitro research use only