

Rat Agrn/Agrin ELISA Kit

Catalog No: #EK19575

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Agrn/Agrin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Agrin
Target Name	Agrin
Other Names	Agrin
Accession No.	P25304
Uniprot	P25304
GeneID	25592;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.98ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Agrn

Note: This product is for in vitro research use only