

Mouse Bmp1/Bone morphogenetic protein 1 ELISA Kit

Catalog No: #EK19576

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Bmp1/Bone morphogenetic protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bone morphogenetic protein 1
Target Name	Bone morphogenetic protein 1
Other Names	BMP-1, Mammalian tolloid protein, mTld, Procollagen C-proteinase, PCP
Accession No.	P98063
Uniprot	P98063
GeneID	12153;

Application Details

Assay Type: Sandwich Detect Range: 0.078-5ng/mL Sensitivity: 0.043ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bmp1

Note: This product is for in vitro research use only