

Bovine CREB1/Cyclic AMP-responsive element-binding protein 1 ELISA Kit



Catalog No: #EK19577

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CREB1/Cyclic AMP-responsive element-binding protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cyclic AMP-responsive element-binding protein 1
Target Name	Cyclic AMP-responsive element-binding protein 1
Other Names	CREB-1,Cyclic AMP-responsive DNA-binding protein,CREB,CREB2
Accession No.	P27925
Uniprot	P27925
GeneID	281713;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:23pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CREB1

Note: This product is for in vitro research use only