

Human CREB1/Cyclic AMP-responsive element-binding protein 1 ELISA Kit

Catalog No: #EK19578

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Human CREB1/Cyclic AMP-responsive element-binding protein 1 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Hu |
| Specificity | Natural and recombinant Human Cyclic AMP-responsive element-binding protein 1 |
| Target Name | Cyclic AMP-responsive element-binding protein 1 |
| Other Names | CREB-1 |
| Accession No. | P16220 |
| Uniprot | P16220 |
| GeneID | 1385; |

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CREB1

Note: This product is for in vitro research use only