

Mouse Adam17/Disintegrin and metalloproteinase domain-containing protein 17 ELISA Kit



Catalog No: #EK19579

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Adam17/Disintegrin and metalloproteinase domain-containing protein 17 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Disintegrin and metalloproteinase domain-containing protein 17
Target Name	Disintegrin and metalloproteinase domain-containing protein 17
Other Names	ADAM 17, TNF-alpha convertase, TNF-alpha-converting enzyme, CD156b, Tace
Accession No.	Q9Z0F8
Uniprot	Q9Z0F8
GeneID	11491;

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 10pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adam17

Note: This product is for in vitro research use only