

Pig ADAM17/Disintegrin and metalloproteinase domain-containing protein 17 ELISA Kit



Catalog No: #EK19580

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig ADAM17/Disintegrin and metalloproteinase domain-containing protein 17 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Disintegrin and metalloproteinase domain-containing protein 17
Target Name	Disintegrin and metalloproteinase domain-containing protein 17
Other Names	ADAM 17, TNF-alpha convertase, TNF-alpha-converting enzyme, CD156b, TACE
Accession No.	O77636
Uniprot	O77636

Application Details

Assay Type: Sandwich Detect Range: 3.12-200ng/mL Sensitivity: 1.56ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ADAM17

Note: This product is for in vitro research use only