

Rat Adam17/Disintegrin and metalloproteinase domain-containing protein 17 ELISA Kit



Catalog No: #EK19581

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Adam17/Disintegrin and metalloproteinase domain-containing protein 17 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Disintegrin and metalloproteinase domain-containing protein 17
Target Name	Disintegrin and metalloproteinase domain-containing protein 17
Other Names	ADAM 17, TNF-alpha convertase, TNF-alpha-converting enzyme, CD156b, Tace
Accession No.	Q9Z1K9
Uniprot	Q9Z1K9
GeneID	57027;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.11ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adam17

Note: This product is for in vitro research use only