

## Human CALD1/Caldesmon ELISA Kit

Catalog No: #EK19582

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CALD1/Caldesmon ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Caldesmon
Target Name	Caldesmon
Other Names	CDM,CAD,CDM
Accession No.	Q05682
Uniprot	Q05682
GeneID	800;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.09ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CALD1

Note: This product is for in vitro research use only