

## Rat Acta2/Actin, aortic smooth muscle ELISA Kit

Catalog No: #EK19583

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Acta2/Actin, aortic smooth muscle ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Actin, aortic smooth muscle
Target Name	Actin, aortic smooth muscle
Other Names	Alpha-actin-2,Actsa,Actvs
Accession No.	P62738
Uniprot	P62738
GeneID	81633;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.172ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Acta2

Note: This product is for in vitro research use only